


Introduction to EEE Design and Project

Coordinator: A/P Yvonne Lam
Email: eyhla@ntu.edu.sg

**Email Contact for administrative matters of uploading files,
makeup arrangements, and other issues:**

**Lau Choon Poh (ECPLAU@ntu.edu.sg) and
Teo Kay Leng, Mark (EKLTEO@ntu.edu.sg)**



Course Structure

- Number of Academic Units - **2 AU**
- Course duration - **36 Hours**
- Course schedule
 - Week **1-2**: Lectures. (**YL/LPH, YL/EG**)
 - Week **3,5**: Design classes – LabVIEW programming and exercise, ELVIS instrumentation and exercise.
 - Week **6**: LabVIEW/ELVIS simulation and **Quiz**.
 - Week **7-12**: Project classes – Project activities including design, development, test and evaluation of subsystems and the integration of hardware and software subsystem modules to form a complete system.
 - Week **13**: Project demonstration and report.



Content ↴

Build Content ↴ Assessments ↴ Tools ↴ Partner Content ↴



[Yvonne Lam](#) ↴

Enabled: Statistics Tracking

This folder contains Lecture notes by Yvonne Lam on Course Information & Electro



[Lee PH](#) ↴

Enabled: Statistics Tracking

This folder contains lecture notes by Prof Lee Peng Hin on Control.



[Erry Gunawan](#) ↴

Enabled: Statistics Tracking

Lecture Notes on all the materials regarding to the lab work and procedures for as



[EE2073 Course Reference](#) ↴



NANYANG
TECHNOLOGICAL
UNIVERSITY
AY2122 S2

Yvonne Lam ↴

Build Content ↴ Assessments ↴ Tools ↴ Partner Content ↴



[Introduction to EE2073](#) ↴

Enabled: Statistics Tracking

Lee PH ↴

Build Content ↴ Assessments ↴ Tools ↴ Partner Content ↴



[Lecture Notes_LPH_EE2073](#) ↴

Enabled: Statistics Tracking

Attached Files: [Lecture Slides on Intro to Control Engineering for EE2073.pdf](#) (895.224 KB)

Erry Gunawan ↴

Build Content ↴ Assessments ↴ Tools ↴ Partner Content ↴



[Saman Notes](#) ↴

Attached Files: [File](#)

This file contains Le

Dr Lam

2 / 64

Reference Materials

EE2073 Course Reference ↴

Build Content ↴ Assessments ↴ Tools ↴ Partner Content ↴



[Assessment Submission Matters](#) ↴

Enabled: Statistics Tracking



[HardwareDatasheets](#) ↴

Enabled: Statistics Tracking



[Lab manual for Week 3-12](#) ↴

Enabled: Statistics Tracking



[Project References](#) ↴

Enabled: Statistics Tracking

Report template; How to upload weekly logbooks; and Turnitin Guide for EE2073 Students

OP275; CA3140; LM380 and THAT_2180_Series_Datasheet

Lab Manual

Data Acquisition Using LabVIEW and ELVIS; Extra Notes on EE2073; [Introduction to LabVIEW](#) and Introduction to NI ELVIS II



NANYANG
TECHNOLOGICAL
UNIVERSITY
AY2122 S2


Dr Lam

3 / 64

Syllabus – Project Classes

(3 Hours per week, week 3 – 13)

- Automatic Volume Control for Audio Amplifier System
- A discrete-time sampled-data digital control system based on LabVIEW and ELVIS



Automatic Volume Control for Audio Amplifier System

Dr Lam

4 / 64

Assessment Scheme

- Design/project practical assessment (**75%** of final mark)
 - Week 3-12: practical skills, understanding, participation, plus ...
 - Week 3-12: results/progress on-site, in weekly eLog uploaded
 - Week 13: project demonstration, Q&A
- LabVIEW/ELVIS assessment (**10%** of final mark)
 - LabVIEW/ELVIS quiz will be conducted in week 6
 - Closed book quiz must be done with the lab computer
 - Students must submit their solutions within the allocated time
- Project report assessment (**15%** of final mark)
 - Individual project report for project activities
 - Content should be consistent with weekly logbook uploaded
 - Subjected to Turnitin originality checking on NTULearn
- Necessary conditions to **PASS** the course:
 - Attendance of at least **80%** in design/project practical sessions
 - Individual project report must be submitted by 5pm Mon of week 14



Dr Lam



5 / 64

eLog system

1. Login to <https://een40003.eee.ntu.edu.sg/elogs>


2. Key in your CITS username and password

(*to access the system from outside NTU, user name to use VPN application as specified by CITS web)




Upload INDIVIDUAL Weekly Logbook

- Before END of each session, upload your logbook file(s) to file server. (Regardless whether all tasks are completed)
- Each file < **4MB** in file size. PDF and WORD are supported.
- Upload all your achievement for the day and no further changes can be made until the next new session.
- Any further changes / revised logbook can only be uploaded at the next new session.
- After upload, download the file to check – ensure your supervisor is able to access.
- **If you fail to upload the logbook by the end of each practical session, you may be deemed to have attained unsatisfactory progress for that particular session.**




Upload INDIVIDUAL Weekly Logbook

- Schematics or snap photos of circuits and instrument setups.
- Screen/window capture of waveforms and run-time outputs.
- Tabulation/graphs of measurement results and theoretical calculations.
- Record of observation key points during experiments.
- LabVIEW/other programs, block diagrams, front panels, etc.



Week 13: Project Report

- Upload softcopy in Word doc to NTULearn for Turnitin originality checking <https://NTULearn.ntu.edu.sg>
- NTULearn can be accessed via NTU Intranet, iNTU (<https://iNTU.ntu.edu.sg>).
- Refer to “Turnitin Guide for EE2073 Students.pdf” for procedure.
- The submission requires each student to attach a similarity report generated by TURNITIN for checking similarity. You need to revise the report if the similarity is more than 20%.
- If requested by your lab supervisor, you may also need to submit hardcopy (double-sided, stapled) to the lab, attached with the first page of Turnitin originality report.
- Submit report (via NTULearn.ntu.edu.sg) by 5pm Mon of week 14.



Lab Locations

- Computer Engineering: S2-B4c-15, or
- Electric Power Research: S2-B5c-02




- Tel: +65 6790 5446
- To spend extra time on the project, arrangement can be made with lab staff in advance – subject to availability.
- There will be no supervision for these arrangements.



Useful Information

- Briefing videos @ youtube
- Lab Sessions @ S2-B4c-15 / S2-B5c-02
- eLog submission: <https://een40003.eee.ntu.edu.sg/elogs>
- Week 6: Quiz (MCQ) – all materials from week 1-5
- Formal Report Submission Deadline: **5pm Mon, Week 14.**
- Formal reports are to be submitted online as Assignment from NTULearn for plagiarism check (**do not submit to eLog system**).
- Makeup should be arranged and performed for the student with any valid reasons, such as absence due to school approved activities or medical reasons.
- The absentee need to contact lab staff for make-up arrangement.
- Absent more than **TWO** lab sessions means **FAIL** of this subject.
- Non-submission of final project report means **FAIL** of this subject.




Download LabVIEW: Please explore!

- Latest download (trial, 45 days validity) available from:
http://www.ni.com/download-labview/?gclid=Cj0KEQjwo7auBRCOtoqn_s-G7aMBEiQAxArNrOzGlQmJoubZtT8FqfouS0IeqCRMmEFNt5PC3zxSzqIaAjy48P8HAQ
- There are many tool boxes that support specific applications.
- There are also many sources for beginner to learn.
- Many working sample programs are also available for each function provided by the software.



Syllabus – Lectures (3 Hours per week, Week 1 – 2)


- Topic 1: Introduction and detailed information on this project.
- Topic 2: Principle of feedback control system (LPH)
- Topic 3: Electronic Circuit and System Design
 - * Electronic circuit and system designs for typical electrical/electronic system functionalities
 - * Using OP AMPS and other commercially available IC building blocks.
- Topic 4: Data Acquisition System and Application (EG/SS)
 - * Data acquisition by sampling for realizing discrete-time signals and sampled-data systems
 - * Application to digital control and signal processing system designs



Syllabus – Design Classes

(3 Hours per week, week 3 – 6)

- LabVIEW as the graphical programming system and application software development environment.
- Students need to learn the basic LabVIEW programming skills from “Introduction to LabVIEW”, which is in the NTULearn folder.
- Data acquisition system, setup and practical applications using ELVIS (Educational Laboratory Virtual Instrumentation Suite).
- Use of LabVIEW and ELVIS for data acquisition application development.



Dr Lam

14 / 64

Project Overview

- Objective: design a controller that automatically adjusts the gain of audio amplifier system to give constant level of audio output.
- Discrete-time sampled-data digital control system design using LabVIEW and ELVIS (Software/Hardware)
- Circuit design and implementation (Hardware)
 - * Voltage controlled amplifier (VCA)
 - * Power amplifier (PA)
 - * Volume unit (VU) meter
- Manual and automatic volume control system design using LabVIEW and ELVIS (Software/Hardware)
- Test, evaluation and demonstration of project using LabVIEW and ELVIS (Software/Hardware)
- Refer [EE2073 Project Manual.pdf](#)





Dr Lam

15 / 64

LabVIEW and ELVIS

- LabVIEW Programming
- Automatic Volume Control



Reminders

- Bring **smartphones** to take photos of your work, which may be used for your weekly logbook.
- Bring **thumb drive** to store your work and the weekly logbook.
- Study the instruction manual before and during class.
- Students are allowed to **take home the breadboard and/ PCB**, which should be carefully kept.
- **Return all boards and components after demo done in week 13.**
- Any problem, please first contact your lab supervisor.
- Safety first! Follow lab rules (proper attire/shoes)

Lab Work in Each Week

- The lab manual provides detailed instructions for students to follow.
- Carefully read the manual and know EXACTLY what to do.
- At the end of instructions for each week, there are open questions to answer, which is OPTIONAL.





Dr Lam

18 / 64

Week 3: LabVIEW Programming

- LabVIEW programs are called Virtual Instruments (VIs) and Each VI has 2 Windows
- Front Panel
 - User Interface (UI)
 - Controls = Inputs
 - Indicators = Outputs
- Block Diagram
 - Graphical Code
 - Data travels on wires from controls through functions to indicators
 - Blocks execute by Dataflow



[Introduction to LabVIEW.pdf](#)



Dr Lam

19 / 64


Week 3: LabVIEW Programming

Loops

While Loops

- * terminal counts iterations
- * Always runs at least once
- * Runs until stop condition is met


While Loop






For Loops

- * terminal counts iterations
- * Run according to input of count terminal

For Loop




Week 3: LabVIEW Exercise



Dr Lam

21 / 64

Week 3: ELVIS Instrumentation+Exercise



Dr Lam


22 / 64

Week 3: ELVIS Instrumentation Hardware Specifications

Introduction to NI ELVIS II.pdf

Oscilloscope

- 16 bit resolution
- 1.25 MS/s single channel, 500kS/s two channel aggregate
- 1 to 1.5 MHz Bandwidth
- 1x and **10x probe**
- ± 10 V input range
- AC/DC coupling
- BNC connection



Digital Multimeter

- Isolated measurements
- 5 digit resolution
- 60 VDC, 20VRms, 2 ADC, 2 Arms, 100MΩ





Dr Lam

23 / 64

Week 3: ELVIS Instrumentation

NI ELVISmx Instrument Launcher

- Digital Multimeter (DMM)
- Oscilloscope (Scope)
- Function Generator (FGEN)
- Variable Power Supplies (VPS)
- Bode Analyzer (Bode)
- Dynamic Signal Analyzer (DSA)



Week 3: ELVIS Exercise

- Amplitude Modulation (AM) – analog modulation scheme where the amplitude of a fixed-frequency carrier signal is continuously varied to represent data

$$y(t) = C \sin(\omega_c t) + M \frac{\cos(\phi - (\omega_m - \omega_c)t)}{2} - M \frac{\cos(\phi + (\omega_m + \omega_c)t)}{2}$$




Figure 1. Time domain of an AM signal






Figure 2. Frequency domain of an AM signal



Week 3: ELVIS Exercise



Week 5: DAQ using LabVIEW/ELVIS



Week 5: DAQ using LabVIEW/ELVIS

Integrated DAQ on ELVIS

Impedance Analyzer


- 0.2 Hz to 35 kHz Range
- NPN, PNP, Diode

Other Analyzers:

- Bode Analyzer
- 2-wire Current Voltage Analyzer
- 3-Wire Current Voltage Analyzer

Prototyping Board

- Updated connections
- Detachable
- User-defined Banana Plugs, BNC, D-Sub connectors



Integrated DAQ

- AI sampling rate 1.25 MS/s single channel, 500kS/s two channel
- 16 bit resolution
- AO 2.8 MS/s update rate
- 24 DIO lines, 15 PFI, 2 CTR


Variable Power Supply

- 10 bit resolution
- 0 to +12V, 0 to -12V
- 500 mA current range


Dr Lam

28 / 64

Week 5: DAQ Exercise: DAQOutputControlVCA




Gauge for VCA: -1.3 – 3 V



Dr Lam

29 / 64

Week 5: DAQ Exercise: DAQInputVUmeter




Week 6: LabVIEW/ELVIS Simulation

- To perform simulation of a volume unit (VU) meter using LabVIEW
- VU meter = envelope detector for audio signals
- The peak envelope detector (capacitor) output V_c for the rectified full-wave input waveform $|V_{in}|$ can be simulated using the simplified dual time constant model as ($\tau_1=0.01$ s and $\tau_2=0.3$ s)

$$V_c(t_{n+1}) = (1 - T_c)|V_{in}(t_n)| + T_c V_c(t_n)$$

$$T_c = \begin{cases} T_{c1} = e^{-\frac{\Delta t}{\tau_1}}, & \text{attack mode: } |V_{in}(t_n)| \geq V_c(t_n) \\ T_{c2} = e^{-\frac{\Delta t}{\tau_2}}, & \text{release mode: } |V_{in}(t_n)| < V_c(t_n) \end{cases}$$

Week 6: LabVIEW/ELVIS Simulation



[VUmeterSim](#)

[VUmeterModel](#)



Dr Lam

32 / 64

Week 6: About the Quiz


- Quiz assessment on LabVIEW programming
- Different MCQs for different students
- **25 minutes to complete 5 questions** by using the computer in the lab



Dr Lam

33 / 64

Week 7: VCA Breadboard



- Gain controlled by V_c (from DAQ O/P):


$$G = e^{-\alpha V_c}$$
- Need dual supplies $\pm 15V$.

Figure 2: VCA Circuit Configuration and Gain Measurement Setup



Breadboard

- Use for your circuits (to bring home)
- The connections of holes are illustrated by the red lines below




ELVIS Prototyping Board

- Use only for supplies, I/O and instrumentation, not your circuits
- The connections of holes are illustrated by the red lines below



Week 7: VCA DAQ Vc




Week 7: VCA Gain Measurement

- Plot Gain (dB) versus V_C
- Deduce V_C sensitivity which is the inverse of the slope

V_{in_pp} (V)	DAQ V_c (V)
(Use matrix no. U1120736A each digit from left to right)	
1. <u>x</u> = 1.1	1.50
1. <u>x</u> = 1.1	1.00
1. <u>x</u> = 1.2	0.50
1. <u>x</u> = 1.0	0.25
1. <u>x</u> = 1.7	0.00
0.1 <u>x</u> = 0.13	-0.25
0.1 <u>x</u> = 0.16	-0.50
0.1 <u>x</u> = 0.11	-1.00

V_{in_pp} (V)	DAQ V_c (V)	V_{out_pp} (V)	$Gain = V_{out_pp}/V_{in_pp}$	Gain (dB)
1	1.5	1.5	1.00	0
1	1	1	1.00	0
1	0.5	0.5	1.00	0
0.3	0.25	0.25	1.00	0
0.3	0	0	1.00	0
0.1	-0.25	-0.25	1.00	0
0.1	-0.5	-0.5	1.00	0
0.1	-1	-1	1.00	0




Week 7: Caution

- Connect up your circuit while ELVIS is OFF.
- Be extremely careful on wire connections, especially for biasing on ELVIS: +15V, -15V, ground, e.g. with only +15V and no/wrong -15V, the IC may get hot quickly.
- Be extremely careful on wire connections, especially for nearby pins.
- Check biasing first if circuit is not working.
- Use consistent wire color if available, e.g. black/blue/brown for ground, yellow/orange for +15V, pink/violet for -15V, red/green/grey/white for DAQ I/O etc.
- Check supply current flow if possible.


Week 8: PA Breadboard

- Provide the required power gain to drive speaker: ~34dB (50)
- Need **single +15V**, pay attention to groundings (7).




[EE2073 Project Manual.pdf](#)

[Read Datasheet LM380.pdf](#)



Week 8: PA final Breadboard




Week 8: PA Freq Response Measurement

- Bode plot: gain and phase response vs freq




Peak Amplitude (Vp)	
(Use matric no. U	1120736A
each different non-zero digit)	
0.02	
0.0 1 = 0.01	
0.0 2 = 0.02	
0.0 7 = 0.07	
0.0 3 = 0.03	
0.0 6 = 0.06	
0.1 1 = 0.11	
0.1 2 = 0.12	
0.1 7 = 0.17	
0.1 3 = 0.13	
0.1 6 = 0.16	




Week 8: Caution

- Connect up your circuit while ELVIS is OFF.
- Be extremely careful on wire connections, especially for biasing on ELVIS: +15V, ground.
- Ground all the ground pins 3,4,5,10,11,12 to help dissipate heat.
- Separate the ground pin 7 from other grounds.
- Switch off ELVIS & PA as soon as measurements done.
- Check biasing first if circuit is not working.
- Use consistent wire color if available, e.g. black/blue/brown for ground, yellow/orange for +15V, pink/violet for -15V, red/green/grey/white for DAQ I/O etc.
- Check supply current flow if possible.



Week 9: VU Meter Breadboard+Soldering

- Based on rectifier concepts.
- Used for measurement of signal amplitude (envelope).
- Capture $VU_{output1}$, $VU_{output2}$ and $VU_{output3}$.



EE2073 Project Manual.pdf

Read Datasheet CA3140.pdf




VU (Volume Unit) Meter

Dr Lam

44 / 64

Week 9: Final circuit



Dr Lam

45 / 64

Week 9: VU Meter Outputs

- Observing input and output waveforms.




Dr Lam

46 / 64

Week 9: VU Meter Response

- Plot Measured Gain versus $V_{U_{\text{input}}}$
- Plot $V_{U_{\text{output}3}}$ versus $V_{U_{\text{input}}}$

VUinput Vin_pp (V)	
(Use metric no. U1120736A each digit from left to right)	
5	
9	= 9.1
4.5	= 4.51
3.5	= 3.52
2.5	= 2.50
1.5	= 1.57
0.5	= 0.53
0.2	= 0.26




Dr Lam

47 / 64


Week 9: Caution

- Connect up your circuit while ELVIS is off.
- Be extremely careful on wire connections, especially for biasing on ELVIS: +15V, ground.
- Check biasing first if circuit is not working.
- Use consistent wire color if available, e.g. black/blue/brown for ground, yellow/orange for +15V, pink/violet for -15V, red/green/grey/white for DAQ I/O etc.
- Check supply current flow if possible.



Week 10: Integration and Testing

- Testing gain controlled by Vc (from DAQ O/P)



Week 10: Integration and Testing


- Plot Gain (dB) versus V_C
- Change frequency and explain

Audio Input (V _{pp})	V_C (V)	Audio Output (V _{pp})	Gain (dB)
1.x	1.50		
1.x	1.00		
1.x	0.50		
1.x	0.25		
0.0x	0.00		
0.1x	-0.25		
0.1x	-0.50		
0.1x	-1.00		

V_{in_pp} (V) (Use matric no. U1120736A each digit from left to right)	DAQ V_C (V)
1. <u>x</u> = 1.1	1.50
1. <u>x</u> = 1.1	1.00
1. <u>x</u> = 1.2	0.50
1. <u>x</u> = 1.0	0.25
0.0 <u>x</u> = 0.07	0.00
0.1 <u>x</u> = 0.13	-0.25
0.1 <u>x</u> = 0.16	-0.50
0.1 <u>x</u> = 0.11	-1.00

Week 10: Final board


- To identify the location of any possible errors, you may need to test each individual sub-system as you did the testings in the previous weeks.



Week 11: Manual Volume Control

- VCA (V_c) = actuator
- VU meter (VUoutput3) = sensor
- No feedback path in the system = open-loop system

EE2073 Project Manual.pdf




Block Diagram Representation of Open-Loop System




Week 11: Manual Volume Control

- Refer to DAQOutputControlVCA and DAQInputVUmeter




Week 11: Manual Volume Control

- Operator-assisted (manual) volume control system
- Adjust volume knob based on VU meter reading (loudness)



Week 12: Automatic Volume Control


- Feedback path thru' controller = closed-loop system



Block Diagram Representation of Closed-Loop System



Week 12: Automatic Volume Control




Step-Up-Down Controller:

$$u(n) = \begin{cases} u(n-1) - \Delta, & e(n) > 0, z(n) < r_0 \\ u(n-1) + \Delta, & e(n) < 0, z(n) > r_0 \end{cases}$$


$$e(n) = r_0 - z(n)$$

where



Δ is the controller parameter variable (step size),
 $u(n)$ is the actuator control voltage for VCA,
 $z(n)$ is the VU meter envelope detector output sensor value
 r_0 is the set point value.



Week 12: Automatic Volume Control




Week 12: Automatic Volume Control



Project Report (Drafting)

- After Week 11, start drafting to submit individual project report.
- Extract from the weekly logbook: VCA, PA, VU meter, integrated system, manual and automatic volume control; schematics or snap photos of final circuits and instrument setups, screen/window capture of final waveforms and run-time outputs, main LabVIEW block diagrams and front panels, etc.
- See slides [7](#) and [8](#) for more information.



Report format



EE2073 Report Template.docx

EE2073 Project Report

Automatic Volume Control for Audio Amplifier System

Student Name: XXX XXX XXX
Matric No.: XXXXXXXXX
Project Group: XXX



School of Electrical and Electronic Engineering
Academic Year 20xx/xx
Semester 2

Dr Lam

60 / 64

Week 13: Project Demo & Report

Formal report due on Monday Week 14 (18 Apr 2022) 5pm



Dr Lam

61 / 64


Week 13: Demo

- The demo should be properly prepared.
- Explain your main achievements, comments on the quality of the work, or any problems, possible further improvements...
- Questions will be asked and both students should provide the answers.



Other References

- More LabVIEW/ELVIS/DAQ, download temp license:
<http://www.ni.com/>
- Basic soldering:
<http://www.aaroncake.net/electronics/solder.htm>
<http://www.youtube.com/watch?v=BLfXXRfRIzY>
- EEE safety video:
http://nocw10.eee.ntu.edu.sg/safety-videos/safety_video_undergraduate.wmv
- EEE courses:
 - EE2001 Circuit Analysis
 - EE2002 Analog Electronics
 - EE2010 Signals and Systems





NANYANG
TECHNOLOGICAL
UNIVERSITY

THE END

